



# MR Hardware and Imperfections

Sebastian Littin

*Medical Physics, University Medical Center Freiburg, Germany*

Ankara, May 26<sup>th</sup>, 2024

# Overview

---

- MRI Hardware Components
- Main Magnetic ( $B_0$ ) Field
- Shims
- $B_1 (+ / -)$  Fields
- Gradients
- System Timing



# MRI Hardware Components

---


What is required for the (N)MR Signal?

What may cause trouble?

- Sample => Signal source ( $^1\text{H}$ ,  $^{13}\text{C}$ ,  $^{17}\text{O}$ , ...)
- Magnet => Polarization (T; static (0Hz); ppm)
- RF coil => Excitation/ Acquisition ( $\mu\text{T}$ ; MHz; xx%)

⇒ Basic NMR Spectrometer

- Shims => Field Homogenization
- Gradient Coils => Signal Localization (mT; mT/m; x%)



# Main Magnetic ( $B_0$ ) Field

- Permanent magnet(s): Low Field
- Electromagnets: Most simple => Solenoid
- Mainly superconductive electromagnets for diagnostic imaging
  - Cryogenic temperature
- Shielding to suppress  $B_0$  outside
- Optimized for ppm (parts per million) accuracy

⇒ Why don't we have perfect  $B_0$ , if the main magnet has ppm accuracy?


(E)



(F)




NMR in Biomedicine, 20 November 2023, DOI: (10.1002/nbm.5052)




# Susceptibility Differences

- Susceptibility ( $\chi$ ): Measure to which extent a substance is magnetized in an external magnetic field




- Substances with different  $\chi$ : Local  $B_0$  gradients
- Cause of  $B_0$  inhomogeneities






# $B_0$ mapping

- Precession frequency depends on  $B_0$ :  $\omega(\vec{x}) = -\gamma B(\vec{x})$
- Phase evolution  $\propto B_0$
- Measure multiple TEs
- Example: Flash sequence

- $B$  large




- $B$  small




# $B_0$ Shimming


- Shimming:  $B_0$  correction/ field homogenization
- Spherical (volume) harmonics
- Resistive coils
- Source of distortion: Inside
- Correction fields: Outside
- Perfect correction not possible
- Temporal stability? Drift?



Handler (2020), Biomed. Phys. Eng.  
Express 6 045022





[https://en.wikipedia.org/wiki/Spherical\\_harmonics](https://en.wikipedia.org/wiki/Spherical_harmonics)



# Radio Frequency Coils

$B_1^+$ : Transmit or Excitation

- Body Coil/ Local transmit coils
- Aim: Homogeneous excitation



Tesfai; JMR (2020) 319,106825


$B_1^-$ : Signal Reception

- Local receiver coil (arrays)
- Localized reception for PI



Siemens Healthineers




⇒ High power vs low signal amplitude



Siemens Healthineers

# RF coils: B<sub>1</sub> Excitation Field Mapping

- B<sub>1</sub><sup>+</sup>: Transmit radiofrequency field
- w/wo Preconditioning RF pulse
- Centric *k*-space reordering  
    ⇒ Residual longitudinal magnetization
- $B_1^+ \propto (S_2/S_1)/\alpha_{prep}$
- Further methods
- $SAR \propto B_0^2 \alpha^2$





# RF coils: Coil Sensitivity Mapping

---

- $B_1^-$ : Localized signal reception
- Important for parallel imaging
- Excitation using Body Coil (BC)
- Multiple receiver Local Coils (LC)

$$\Rightarrow B_1^- \propto \frac{LC}{BC}$$


- Numerous software packages available




# Gradient Coils for Spatial Localization

---


- Linearly varying magnetic field ( $B_z$ )  
= constant gradient
- Local variation of the Larmor frequency
- Linear relation between object and frequency space
- Straight-forward Fourier Transform (FFT)
- Usually three orthogonal magnetic field gradients along spatial directions x,y and z




# How linear are our encoding fields?



clinical scanner, XZ (coronal) slice  
FOV=50cm, iso-center



FOV=50cm, iso-center  
*no gradient distortion correction*




FOV=50cm, 20cm above iso-center  
*no gradient distortion correction*

⇒ No gradient system is perfect

# Gradient Coils in practice

---


- High gradient amplitude (mT/m, scales with current)
  - Fast switching/ slew rate (T/s, scales with voltage)
  - Modern PWM amplifiers >1000A; >2kV
  - Currents inside a magnetic field:  $F = I \int d\mathbf{l} \times \mathbf{B}$ 
    - Similar to a loudspeaker: Source of most noise in MRI
    - Multiple kW of heat: active cooling required
  - Fast switching:
    - „Eddy currents“: Opposing fields from currents induced in cryostat with decay times > 1s
- ⇒ Magnetically screened, „shielded“ coils
- In practice: Limitations given by physiology (PNS)






# Gradient Mapping

---

- Distortion phantoms
- Dedicated devices  
(Hall probes)
- Field probes



ISMRM 2022 #5032: Distortion Measurements Using a Large-volume, Lightweight, Modular MRI Phantom With Solid Signal Sources








<https://www.sunnuclear.com/products/large-field-mr-distortion>






# Gradient Mapping of Insert Coils

- Play out blips on additional gradient channel(s)
- GRE sequence to measure phase w/wo blips
- Calculate phase difference



# Eddy-current mapping

- Switch on gradient, hold and wait (minimize effects from switching on)
- Switch-off fast and measure
- Map phase evolution from eddy-currents (phase differences similar to field mapping)



# Before you start on the scanner



---

- Peripheral nerve analysis
- Avoiding acoustic resonances (for fast sequences)



# Gradient Interactions: Safety I (PNS)

- Unwanted Peripheral Nerve Stimulation (PNS)
- Maxwell-Faraday equation:  $\nabla \times E = -\frac{\partial B}{\partial t}$
- Voltage along nerve fibers, strong enough to trigger discharge through ion channels
- Usually below cardiac stimulation thresholds
- Magnetic fields generated by gradient coils go beyond FOV
- All field components ( $B_x, B_y, B_z$ ) interact with the human body
- Prediction in Pulseq possible



[1] M Davids et. Al. Scientific Reports 2017;7:5316.

[2] M Davids et Al. MRM 2019;81:686-701.

# *Pulseq* PNS prediction

---

- Peripheral neural stimulation (PNS) prediction on Siemens is based on the so-called SAFE model  
F.X. Herbank and M. Gebhardt. SAFE-Model - A New Method for Predicting Peripheral Nerve Stimulations in MRI. ISMRM 2000, #2007.  
<https://cds.ismrm.org/ismrm-2000/PDF7/2007.PDF>
- Open-source implementation by Filip Szczepankiewicz and Thomas Witzel:  
[https://github.com/filip-szczepankiewicz/safe\\_pns\\_prediction](https://github.com/filip-szczepankiewicz/safe_pns_prediction)
- Direct interface in Matlab-Pulseq



# Caveat & Solution


---

- Need scanner-specific model parameters
- Extract your gradient system description file from IDEA
  - Type “sys” in the IDEA shell selecting your system  
(you may need to select something else first to see the verbose output)
  - Note the strings after “GPA Type” and “GC Type”
  - Go to `C:\MIDEA\N4_VE####\n4\pkg\MrServers\MrMeasSrv\Config\InitMeas` and pick the file named `MP_GPA_<your_GPA>_<your_coil>.asc` and copy it somewhere where your Matlab or Python can read it
  - *Path under NumarsX is left as an exercise for the interested reader*
- You can now predict PNS in the same way IDEA and scanner do it

# Example: DW-EPI

---


- bFactor=1000, TE=78ms, Gmax=38mT/m SR=180T/m/s



# PNS Stimulation: DW-EPI




---

```
seq.calcPNS('idea/asc/MP_GPA_K2309_2250V_951A_AS82.asc');
```




# Gradient Interactions: Safety II (Mechanical)


- Mechanical:  $F = I \int d\mathbf{l} \times \mathbf{B}$ 
  - Vibrations
  - Mechanical resonance
  - ⇒ Acoustic noise
  - Electromechanical coupling
  - ⇒ Helium boil-off
- Prediction of “Forbidden Frequencies” possible  
(additional knowledge about hardware required)



R Kimmlingen. MAGNETOM Flash (68) 2/2017



SA Winkler et. Al. NeuroImage 2018;168:59 –70.



# Acoustic resonance analysis


---

- Extract your gradient system description file from IDEA as above
- Create your sequence in memory  
(populate the “seq” object with events)
- Run “gradSpectrum.m” script in “matlab/demoUnsorted”
  - Set ascName=‘....’ to your gradient system .asc file to see resonances marked



# Acoustic analysis: DW-EPI

---




seq.sound(); % another option for acoustic assessment in Matlab ;-)



# System Timing

---

- Different sub-systems have to operate synchronously
  - RF / ADC (MHz) / Gradients (kHz)
  - Different orders of magnitude  
⇒ RF pulse synced with gradient?
  - ⇒ ADC center at TE?
- Nominal vs real trajectory
- Measurements, eg field probes  
⇒ Skope uses Pulseq...



# Acknowledgements

*University Medical Center Freiburg*

Whole Pulseq Team

All current and previous  
colleagues and collaboration  
partners



Funded by  
the European Union

HORIZON-WIDERA-2021-ACCESS-03  
MRITwins

